



## 4 Easy Steps!




(1) Collect Specimen




(2) Remove Cap & Dip Test Strips for 30 Seconds

(3) Return the Cap & Wait 3 Minutes



(4) Read Test Lines for Results

(one line = **negative**, two lines = **positive**)



*This test provides a presumptive diagnosis for pregnancy. A confirmed pregnancy diagnosis should only be made by a physician after all clinical laboratory findings have been made. See package insert for complete details.*

**QuikTest USA, Inc. – 10959 Handel Place Boca Raton, FL 33498 – (866) 989-9300 – [www.quiktestusa.com](http://www.quiktestusa.com)**